



[www.Duramedic.eu](http://www.Duramedic.eu)

# The 'Protect-a-Med' Disposable Linen Range from Duramedic™



**Critical Healthcare Ltd, Comagh Business Park  
Kilbeggan, Co. Westmeath, Ireland**

Tel: +353 (0) 57 93 33111 Fax: +353 (0) 57 93 33115

**These products have passed FDA inspection, ISO 9001:2000, ISO 13485:2003, and CE.**


**Our products have also earned FDA 510(k) approval.**



[www.Duramedic.eu](http://www.Duramedic.eu)

## BL3 I0 Protect-a-Med Disposable Linen Day Pack

(R.R.P. €85/Box of 20 packs)



- An effective barrier to bacteria
- Drives down or eliminates the cost of cleaning and or replacing stretcher mattresses
- Significantly reduces the risk of cross infection to staff and patients
- Reduces the time spent cleaning the vehicles, thus minimizing down time of the ambulance and staff
- Reduces risk of staff contamination from packing and transporting soiled laundry.



Supplied in a perforated polythene bag: 210(w) x 260(h) x 40(d) mm. Packed weight approx. 220g.

The day pack is sold and packaged in cartons of 25 packs.

### **Contents:**

- I x BLI40 Disposable Pillow Case
- I x BLI15 Disposable Fitted Stretcher Sheet
- I x BLI20 Disposable Patient Blanket




[www.Duramedic.eu](http://www.Duramedic.eu)

## **BLI40 Protect-a-Med Disposable Pillow Case**

(R.R.P. €64/Box of 100)



- An effective barrier to bacteria
- Drives down or eliminates the cost of cleaning and or replacing stretcher mattresses
- Significantly reduces the risk of cross infection to staff and patients
- Reduces the time spent cleaning the vehicles, thus minimizing down time of the ambulance and staff
- Reduces risk of staff contamination from packing and transporting soiled laundry.



- Soft Fluid Resistant Linen
- Strong Durable and Lightweight
- Non-Woven and Breathable
- Colour: White
- Dimensions 530(w) x 750(h) mm. Weighs approx. 37g




[www.Duramedic.eu](http://www.Duramedic.eu)

## **BLI I5 Protect-a-Med Disposable Fitted Stretcher Sheet**

(R.R.P. €65/Box of 50)



- An effective barrier to bacteria
- Drives down or eliminates the cost of cleaning and or replacing stretcher mattresses
- Significantly reduces the risk of cross infection to staff and patients
- Reduces the time spent cleaning the vehicles, thus minimizing down time of the ambulance and staff
- Reduces risk of staff contamination from packing and transporting soiled laundry.



- Soft Fluid Resistant Linen
- Strong, Durable and Lightweight
- Non-Woven and Breathable
- Ideal for use as a transfer sheet
- Colour: Light Blue
- Dimensions 750(w) x 1800(h) mm. Weighs approx. 96g




[www.Duramedic.eu](http://www.Duramedic.eu)

## **BL210 Protect-a-Med Disposable Patient Blanket**

(R.R.P. €50/Box of 50)



- An effective barrier to bacteria
- Drives down or eliminates the cost of cleaning and or replacing stretcher mattresses
- Significantly reduces the risk of cross infection to staff and patients
- Reduces the time spent cleaning the vehicles, thus minimizing down time of the ambulance and staff
- Reduces risk of staff contamination from packing and transporting soiled laundry.



- Fluid Resistant Linen
- Strong, Durable and Lightweight
- Non-Woven and Breathable
- Ideal for use as a privacy sheet
- Colour: Dark Blue
- Dimensions 1000(w) x 2100(h) mm. Weighs approx 81g



[www.Duramedic.eu](http://www.Duramedic.eu)

## **BL910 Protect-a-Med Fleece Rescue Blanket**

(R.R.P. €6.95)



- The ideal replacement for existing expensive ambulance blankets
- A Warm, Waterproof Protective Fleece which is Strong and Durable yet Lightweight
- The Fleece Rescue Blankets strength makes it Ideal for use as a carry sheet
- Cost effective - can be used as a single use item



- Dimensions 2290(h) x 1520(w) mm.
- Weighs approx 690g (packed)
- Colour: Dark Blue

### **BL910 Fleece Rescue Blanket - Independent Laundering Trial**

The Central Laundry in St. Loman's Hospital have independently tested the BL910 Fleece Rescue Blanket and have confirmed that in their opinion it is suitable for use as a reusable/laundryable Ambulance Blanket. They washed and tumble dried the Fleece Rescue Blanket 20 times at normal-high temperatures (2 washes at 40°C, 14 washes at 80°C and 4 Washes at 90°C) with no significant deterioration in the properties of the blanket. The Central Laundries opinion was that the blanket would withstand 40 such washes before deterioration.